

# Gamma Camera SPECT Uniformity

 Project

8th January - 22nd April, 2024

Word count: 3800 (excl. appendices)

**Abstract:** Gamma cameras are imaging tools used in hospitals to view physiological processes by imaging the gamma radiation emitted by radio-pharmaceuticals. As a part of quality control measures to ensure the functionality of gamma cameras, uniformity tests need to be completed to ascertain whether the devices are functioning correctly. This involves imaging a uniform body of radiation, also called a test phantom, and checking that the DICOM images produced by a gamma camera show a similarly uniform image. Traditionally, medical professionals have inspected DICOM images and measured the uniformity across limited parts of the image. We have aimed to create a piece of software as a proof of concept that allows uniformity to be calculated across an entire DICOM image at once, that has the potential to increase the efficiency and reliability of quality control tests for gamma cameras.

Lyall Stewart

Leah Wells

Zac Baker

Alex Toogood-Johnson

Josh Scates

 School



## Contents

<b>1</b>	<b>Introduction and Background</b>	<b>3</b>
1.1	Introduction to Nuclear Medicine . . . . .	3
1.1.1	What is Nuclear Medicine? . . . . .	3
1.1.2	What is a Gamma Camera? . . . . .	3
1.1.3	How Does a Gamma Camera Work? . . . . .	4
1.2	SPECT . . . . .	4
1.3	Quality Control . . . . .	5
1.4	Uniformity . . . . .	5
<b>2</b>	<b>Methods</b>	<b>7</b>
2.1	DICOM files . . . . .	7
2.1.1	The DICOM File Format . . . . .	7
2.1.2	Obtaining Metadata . . . . .	8
2.1.3	Obtaining Image Data . . . . .	9
2.2	What is a Convolution? . . . . .	10
2.3	Development Tools . . . . .	10
2.3.1	Programming language . . . . .	10
2.3.2	GUI framework . . . . .	11
<b>3</b>	<b>Results</b>	<b>12</b>
3.1	Final Program . . . . .	12
3.2	Testing . . . . .	14
<b>4</b>	<b>Discussions</b>	<b>15</b>
4.1	Execution Speed . . . . .	15
4.2	Code Structure Improvements . . . . .	15
<b>5</b>	<b>Acknowledgements</b>	<b>16</b>
<b>6</b>	<b>References</b>	<b>16</b>
<b>A</b>	<b>Code Appendix</b>	<b>18</b>
A.1	dicom_functions.py . . . . .	18
A.2	uniformity_functions.py . . . . .	21
A.3	gui.py . . . . .	23
A.4	config.json . . . . .	34
A.5	dicom_elements.json . . . . .	35

## 1 Introduction and Background

The RD&E Nuclear Medicine Department currently performs various types of quality control tests on Gamma cameras [3]. One of these tests involves calculating a measure of uniformity over a tomographically acquired SPECT image. Our task was to write a piece of software that allows a user to find the differential and integral uniformity of a DICOM image, by accepting an image from the user, applying a 9-point weighted convolution filter to smooth out the pixels, and using the formulae specified in the NEMA handbook provided to us to obtain the uniformity percentage values.

### 1.1 Introduction to Nuclear Medicine

#### 1.1.1 What is Nuclear Medicine?

Nuclear medicine is a specialised field of medicine that uses radioactive tracers known as radio-pharmaceuticals to assess bodily functions and to diagnose and treat disease.” Approved tracers are called radio-pharmaceuticals since they must meet the FDA’s exacting standards for safety and appropriate performance for the approved clinical use” [3]. Some commonly used tracers used in nuclear medicine are Technetium-99m, chosen due to its short half-life of 6 hours, and Iodine-131 [2].

#### 1.1.2 What is a Gamma Camera?

Gamma cameras are imaging tools used in hospitals to view physiological processes by imaging gamma radiation emitted by radio-pharmaceuticals [4]. A gamma camera produces SPECT images (see Section 1.2) from the gamma ray emissions, and these can be analysed to diagnose a condition. For example, a high radiation count in a specific area of the body, when the radio-pharmaceuticals have been attached to substances such as octreotide or dotatate, can be indicative of the presence of neuroendocrine tumours [15].



Figure 1: A modern dual-head gamma camera [9]

### 1.1.3 How Does a Gamma Camera Work?

A radio-pharmaceutical such as Technetium-99m can be attached to a drug so that it will target a specific area of the body that needs to be imaged [2]. As mentioned in section 1.1.2, drugs such as octreotide and dotatate can be attached to these radio-pharmaceuticals so that the radioactive sources are drawn to tumour cells, although many other possibilities exist [15].




Figure 2: Gamma camera diagram [16]

Gamma cameras detect high-energy photons released from the radioactive atoms in a patient's body. These photons escape from the patient's body from all directions and are passed through a collimator. This ensures that only photons travelling vertically upwards will be detected, which reduces noise and measures with more accuracy the density of the radiation emitted from different areas. The photons then hit a scintillation crystal, which converts them into light. This light then has to be amplified in a photo-multiplier tube array and processed by a computer to produce a 2D image of the distribution of photons. The camera takes these images from 360°, so they can be layered to form a 3D image of the organ or process you are interested in [4]. The areas that you are interested in will have a higher concentration of photons.

## 1.2 SPECT

A SPECT scan is a type of imaging that uses a radioactive substance and a detector to produce 3D images to show physiological processes. A gamma camera produces SPECT scans [9]. When multiple 2D projections are acquired of an object from multiple angles, mathematical tools can then be used to reconstruct a 3D image of that object.

When using a gamma camera, the image reconstruction is created from planar (2D) images from every angle, which shows the number of gamma rays detected [10].

### 1.3 Quality Control

Quality control means ensuring that the level of performance of a system remains at a stable operating level within the most recent performance limits. A good quality control system ensures that the people closest to the work have the knowledge, tools, and procedures in place to review, reflect, and react, with immediacy, to discrepancies between observed and agreed levels of performance [6]. This is necessary because an error in the system could lead to the misdiagnosis of patients. In the case of gamma cameras, it could also lead to increased exposure to radiation if the test must be redone.

### 1.4 Uniformity

The uniformity of a gamma camera is normally defined as the difference between the maximum and minimum counts per field element (or count density) expressed as a percentage of the mean number of counts per element or of the average of the maximum and minimum values [8]. For example, the uniformity of a gamma camera is often tested using a cylinder of radiation, if the gamma camera is working perfectly this will produce one section of equal maximum density surrounded by equal minimum density; in reality this is not the case, so uniformity tests must be used to confirm that the variance in count density is due to background noise and not a fault with the gamma camera. There are two types of uniformity: differential uniformity and integral uniformity. If the value of differential uniformity exceeds 4% on a daily test, then it should be retested and considered as potentially being faulty [12].

$$\text{Integral uniformity} = \pm 100 \times \frac{(max - min)}{(max + min)} \quad (1)$$

From the NEMA handbook given to us for this project: "For pixels within each area (CFOV and UFOV), the maximum and minimum values are identified from the smoothed data. The difference between the maximum and the minimum is divided by the sum of these two values, and multiplied by 100." Here, CFOV and UFOV stand for Central Field Of View and Useful Field Of View. As we have cropped each DICOM image to a specific width, length, and height, we do not need to regard the fields of view, and can instead focus on values for the whole image. Thus, to find the integral uniformity of a DICOM image, we just smooth it with a convolution filter, crop it to the preferred size, and substitute the maximum and minimum pixel values across the whole image into the formula above.

$$\text{Differential uniformity} = \pm 100 \times \frac{(max - min)}{(max + min)} \quad (2)$$

From looking at the equation, it may seem that differential uniformity is the same as integral uniformity, however, the main difference here is which values are chosen as *min* and *max* to substitute into the equation. From the NEMA handbook: "The smoothed data are treated as a number of rows (X slices) and columns (Y slices) Each slice is processed by starting at the beginning pixel for the respective FOV. A set of 5 contiguous pixels is examined to find the minimum and maximum pixels. The differential uniformity is calculated using these values. The next set of 5 pixels is analysed by stepping forward one pixel and again determining the percent uniformity. This is repeated until the outermost pixel is reached. The maximum differential uniformity is found in the slice. This process is then repeated for all the slices, and the maximum value for all slices is reported." To summarise, each DICOM image is made of slices, e.g. a cropped, smoothed 8x8x8 DICOM image has 8 layers, each of which are 8 pixels wide by 8 pixels long.

7	8	4	9	6	3	10	9
5	4	6	2	3	7	6	7
5	7	3	1	5	4	2	3
3	2	9	9	6	7	8	6
2	3	4	3	10	11	9	6
7	2	3	5	6	2	3	4
1	4	1	2	7	5	4	5
5	4	6	2	3	7	6	7

For each of those slices in the image, take a set of 5 contiguous values, e.g.:

7	8	4	9	6
---	---	---	---	---

We can find the values of  $max$  and  $min$  from these values, being 9 and 4, and substitute these into the equation for differential uniformity to get a uniformity value of 38.46%, to 2 decimal places. Let us step forward one pixel and consider the values:

8	4	9	6	3
---	---	---	---	---

Here, we can take the values for  $max$  and  $min$  as 9 and 3, and get a differential uniformity of 50%. We would then repeat this for every row and column in the layer, 5 pixels at a time, to get a uniformity value for that layer (this being the maximum value for any 5 pixels in that layer). The maximum value out of all of the layers in the image is taken as the differential uniformity for that image. Due to each calculation requiring 4 additional pixels to the right or down, this leaves a 5 pixel thick band around the left and bottom of the circle in each layer, which cannot be calculated with and are ignored.

## 2 Methods

### 2.1 DICOM files

At first, we decided to use the Python package *pydicom* in order to read DICOM files, however, after that we decided to instead research DICOM files and the DICOM specification to try to manually parse the binary of a DICOM file in order to get the metadata and image data. We eventually succeeded in doing this to the most extent, although there are still a few flaws in the code.

#### 2.1.1 The DICOM File Format

The results of gamma camera scans, as standard in medical imaging, are stored in the Digital Communications and Information in Medicine (DICOM) file format, allowing cross-compatibility with a number of different software tools. While there are many open-source, existing tools to parse these, which we did use at the start of the project, we later decided to do this step ourselves.

We researched the structure of DICOM files and found out that there are 2 main parts: the header and the dataset. The header is made of a preamble and a prefix. Firstly, the preamble is a sequence of 128 bytes, which are usually all set to 'OOH'. The prefix is a sequence of 4 ASCII-encoded bytes, which stand for the letters 'D', 'I', 'C' and 'M'. Thus, we ignore the first 132 bytes of any DICOM file [10].




Figure 3: Structure of a DICOM file [10]

The dataset is comprised of repeating data elements, which are split further into either

metadata or image data elements. Each element has 5 parts : the group ID, the element ID, the VR (Value Representation), the group length, and the element data. The group ID, the element ID, and the VR are always 2 bytes each. The size of the value length field is determined by the which VR is used, and the value of the element data field is determined by the data length field. There are 34 different VRs, each of which are a 2-byte ASCII-encoded string. Each VR holds a different type of data, for example the VR 'UL' is an unsigned 32-bit binary integer. The full list of VRs can be found in the DICOM specification [6].

Group Number	Element Number	Value Representation	Value Length	Value Field
2 bytes	2 bytes	2 bytes	2 bytes	"Value Length" bytes

Figure 4: Structure of a DICOM element [11]

Each group ID and element ID pairing, which take the form of 2 4-digit numbers e.g. (0100,0010), correspond to a specific tag, the full list of which, can be found from pages 23 to 183 of the registry of DICOM data elements [8]. In the image above, the VR is neither OB, OW, SQ or UN, so the value length is 2 bytes, however, if the VR is one of those 4, then a 4 byte value length shall be preceded by 2 reserved bytes, which can be ignored when we parse the file, meaning 6 bytes are used instead of 2. One of the challenges we faced when manually decoding the binary was when we encountered one of these 4 VRs. SQ, which stands for Sequence, has multiple data elements contained within, and we found this very difficult to parse.

After the metadata elements, there is the image data element. This stores the pixel data for the entire DICOM image, and we found in the DICOM files we used that this had a tag beginning '7FEO' [7]. We initially had trouble decoding this, but upon looking more closely at the DICOM specification, we eventually figured out how to convert the binary into a Python *numpy.ndarray* representing the image.

### 2.1.2 Obtaining Metadata

Our main success with manually parsing the DICOM files was obtaining most of the metadata for a DICOM file, as we could get the group ID, element ID, VR, length and data field for almost all of the metadata elements. We also web-scraped the registry of DICOM data elements into a JSON file, so that we could convert the group-element pair into a tag [5]. Comparing our result with the metadata output of the *pydicom* package, we found out we were almost completely accurate, with our only issue being that our code couldn't decode any SQ (Sequence) elements.



```

File Meta Information Group Length 02000000 212
File Meta Information Version 02000100 b'\x00\x01'
Media Storage SOP Class UID 02000200 1.2.840.10008.5.1.4.1.1.20
Media Storage SOP Instance UID 02000300 1.2.840.113619.2.184.31108.192168116
Transfer Syntax UID 02001000 1.2.840.10008.1.2.1
Implementation Class UID 02001200 1.2.840.113619.6.281
Implementation Version Name 02001300 Xeleris 3.1108
Source Application Entity Title 02001600 XELERIS-R
Group Length 00080000 544
Image Type 00080008 DERIVED\PRIMARY\RECON TOMO\EMISSION
Instance Creation Date 00080012 20231124
Instance Creation Time 00080013 152143.0000
Instance Creator UID 00080014 1.2.840.113619.6.281
SOP Class UID 00080016 1.2.840.10008.5.1.4.1.1.20
SOP Instance UID 00080018 1.2.840.113619.2.184.31108.192168116
Study Date 00080020 20230728
Series Date 00080021 20231124
Acquisition Date 00080022 20230728
Content Date 00080023 20231124
Study Time 00080030 103856.00
Series Time 00080031 152054.00
Acquisition Time 00080032 103856.00
Content Time 00080033 152054.00
Modality 00080060 OT
Manufacturer 00080070 GE MEDICAL SYSTEMS
Institution Name 00080080 Royal Devon and Exeter Hospital
Station Name 00081010 XELERIS-R
Study Description 00081030 User&Basic&Tomo
Series Description 0008103e Volumetrix MI RESULTS AC
Manufacturer's Model Name 00081090 INFINIA
Private Creator 00090000 470
Private Tag Data 00090010 GEMS_GENIE_1
UNKNOWN 00091001 GEMS_GENIE
Study Name 00091010 Tomo
Study Flags 00091011 17825792
Study Type 00091012 0
Dataset UID 0009101e 1.2.840.113619.2.184.31108.192168116
Series Object Name 00091020 Tomo
Series Flags 00091021 16777216
Initiation Type 00091023 0

```

Figure 5: Part of the metadata obtained through manually parsing the binary

```

Group Length: 4194316
Image Type: ['DERIVED', 'PRIMARY', 'RECON TOMO', 'EMISSION']
Instance Creation Date: 20231124
Instance Creation Time: 152143.0000
Instance Creator UID: 1.2.840.113619.6.281
SOP Class UID: 1.2.840.10008.5.1.4.1.1.20
SOP Instance UID: 1.2.840.113619.2.184.31108.192168116118.1700839067.58997240
Study Date: 20230728
Series Date: 20231124
Acquisition Date: 20230728
Content Date: 20231124
Study Time: 103856.00
Series Time: 152054.00
Acquisition Time: 103856.00
Content Time: 152054.00
Accession Number:
Modality: OT
Manufacturer: GE MEDICAL SYSTEMS
Institution Name: Royal Devon and Exeter Hospital
Referring Physician's Name:
Station Name: XELERIS-R
Study Description: User&Basic&Tomo
Series Description: Volumetrix MI RESULTS AC
Name of Physician(s) Reading Study:
Operators' Name:
Manufacturer's Model Name: INFINIA
Private Creator: 1052
Private tag data: GEMS_GENIE_1
[Unknown]: [(0013, 0010) Private Creator                                LO: 'GEMS_GENIE_1'
(0013, 1014) [Original Image Num]                                     SI: Array of 128 elements]
[Study Name]: Tomo
[Study Flags]: 4097
[Study Type]: 0
[Dataset UID]: 1.2.840.113619.2.184.31108.192168116118.1700839067.58997240
[Series Object Name]: Tomo
[Series Flags]: 1
[User Orientation]:
[Initiation Type]: 0
[Initiation Delay]: 0
[Initiation Count Rate]: 0

```

Figure 6: Part of the metadata obtained via the *pydicom* package

Here in the images above and below, you can see a list of the metadata that you would find in a standard DICOM file. There are elements such as the acquisition date, the manufacturer of the gamma camera, and the name of the institution who owns the gamma camera. We had just over 5,200 entries in our JSON file, which meant our program was capable of identifying that many elements of metadata, however there are still a few unknown group-element ID pairings, which have just been left as 'UNKNOWN' in the metadata output.

Although our program did not end up using a majority of the metadata we obtained through reading the binary, we feel that coding this helped us gain a better understanding of the structure of DICOM files, and it helped us understand how to get the image data from the binary, which our program did use.

### 2.1.3 Obtaining Image Data

Once we had managed to obtain the metadata from the SPECT DICOM files that we were given, we found it easy to gain the raw data for the image, that being a Python list some 4,194,304 bytes in length, however it then took us quite a long time to figure out how to convert this into an image, although in the end it turned out to be relatively easy.

The metadata with group (0028, 0100) contains the value for the bits allocated per pixel [7]. For all the DICOM files we tested, this turned out to be 16 bits, or 2 bytes per pixel. From here, the groups with ID (0028, 0008), (0028, 0010) and (0028, 0011) in the metadata contained the values for the number of rows, columns and planes in the SPECT image. From this point, we converted the Python list to a *Numpy* array, and then used the *Numpy* *.reshape()* function to convert the *Numpy* array into a 3D array with the same dimensions as the image.

## 2.2 What is a Convolution?

"A convolution is essentially a filter that can be applied to a set of pixels in order to emphasize or smooth a particular feature. A convolution uses a matrix called a kernel, and each pixel in the resultant image is a function of the nearby pixels, and the input image, with a kernel as the function" [5]. Mathematically :

$$g(x, y) = \omega(x, y) * f(x, y) = \sum_{i=-\infty}^{\infty} \sum_{j=-\infty}^{\infty} \omega(i, j) f(x - i, y - i)$$

Where  $g(x, y)$  is the filtered image,  $f(x, y)$  is the original image, and  $\omega$  is the filter kernel [14]. For more information, we recommend watching the 3Blue1Brown YouTube video on convolutions [1]. Different kernels can create different effects, such as sharpening or blurring an image, but the kernel we were given is intended to smooth a layer of a DICOM image. In the NEMA handbook provided to us, we were given the following kernel to use.

1	2	1
2	4	2
1	2	1

In order to understand how a convolution kernel works, let us take a small section of the arbitrary layer in section 1.4, focusing specifically on this '6' pixel:

7	8	4	9
5	4	<b>6</b>	2
5	7	3	1
3	2	9	9

In order to obtain the new value for the '6' pixel, we must first multiple each pixel in the 3x3 section centered around '6' by the corresponding pixel in the convolution kernel, obtaining the following matrix:

8	8	9
8	24	4
7	6	1

From this, we can take the mean value of this matrix, 8.33, to 2 significant figures. This is the new value for the '6' pixel. This process is repeated for every pixel in the layer, and every layer in the DICOM image, to get the smoothed image. We then crop the image to a cylinder, discarding the outer regions.

## 2.3 Development Tools

### 2.3.1 Programming language

We decided to use Python in order to code our program, as it is a language that we were all familiar with before starting the project. We did contemplate using a combination of JavaScript, HTML and CSS for a web based design, but we decided the extra work of learning components of a new programming language in addition to completing the project would have been too much to complete. Python also came with the most pre-written libraries involving DICOM files, however, we decided not to use these.

### 2.3.2 GUI framework

Some of the GUI frameworks we considered using were.

- Tkinter
- Kivy
- PyQT5
- CustomTkinter [13]

In the end we decided that Kivy and PyQT5 had a steep learning curve, as we did not previously know how to use them, which left us with a choice between Tkinter and CustomTkinter. Although CustomTkinter had slightly different functionality, it was sufficiently similar to Tkinter, which we had used before, that we felt confident in using it. The advantages of CustomTkinter over Tkinter are a nicer-looking GUI, support for light / dark mode and colour themes, and dynamic scaling when resizing an app window.

### 3 Results

<https://github.com/AlexToogood-Johnson/Gamma-Camera-Uniformity>

This is the GitHub repository where we have put all our code. If you have Git installed on your system, you can run the command:

```
git clone github.com/AlexToogood-Johnson/Gamma-Camera-Uniformity.git
```

Which will download all the files, or else you can download them manually. The appropriate Python libraries can be installed with `pip install -r requirements.txt`.

#### 3.1 Final Program




Figure 7: Our program showing a smoothed and cropped DICOM image

Note the following features:

- **Left & Right:** Increases or decreases the index of the layer shown in the preview
- **Open & Close:** Select a DICOM file for processing, or reset the program to the beginning ready for a new file.
- **Save Layer:** Export the current layer as shown in the preview to an image file.
- **Smooth:** Applies the convolution to the image, as specified in Section 2.2.
- **Reduce:** Removes the top and bottom  $n$  layers of the image, by default will remove all but the middle 40.
- **Revert:** Undo all processing to the DICOM file.




Figure 8: Our program showing the result of a uniformity calculation




Figure 9: Our program's settings page

Again note the following:

- **Crop Size:** affects the number of layers that will be removed from the top and the bottom.
- **FoV Radius:** Controls the radius of the circle which will be considered for uniformity.
- **Step:** See Section 4.1. A step of 1 calculates the uniformity for every layer, a step of 2 considers every other layer etc. Higher values will execute faster, but will be less accurate.

### 3.2 Testing

We were provided with two sample DICOM scans of standard phantoms for testing, from which we obtained the following results:

File Name	Differential Uniformity (%)	Integral Uniformity (%)
TOMO_IRNC001_DS.dcm	3.37	30.14
TOMO_IRChangAC001_DS.dcm	2.70	23.92

Table 1: Testing values for provided DICOM images.

The values for differential uniformity are within the acceptable  $< 3 - 4\%$  range (although, the first value is the high level of what would be acceptable), however, calculated values for integral uniformity seem particularly high, with a large range [12]. In absence of an online source confirming what an acceptable threshold for integral uniformity would be, we have decided to accept the calculated values as accurate. Existing sources only calculated the integral uniformity over small regions of interest, so values they suggest as maximum thresholds are not applicable to us. Because we consider the whole image, our values are expected to be much larger.

## 4 Discussions

### 4.1 Execution Speed

Uniformity calculations, particularly for differential uniformity, are highly computationally intensive tasks. This, without proper optimisation, can result in very slow execution times, on average approximately 2.5 minutes per DICOM image.

This was unacceptably slow, and at 3.65 seconds per layer (for the central 40 layers considered by the program), this is likely slower than manually inspecting each layer for major anomalies. To optimise this, we adjusted the code to allow for multiprocessing, where different calculations could be executed simultaneously by taking advantage of modern multi-core CPUs. This immediately gave an approximately 85% reduction in execution speeds with no impact on the final calculated uniformity values.

To further optimise the program, we decided to only calculate the uniformity for every other layer in the cropped region. Our reasoning was that any significant anomalies would likely affect multiple layers and therefore only considering intermediate layers would be sufficient. This gave us a further approximately 39% reduction when tested. While this did not affect the final uniformity values for the two tested images, this was by chance and it would cause value changes for some images. However, consecutive layers tend to have similar values, so it was decided that this loss of information was justified by the improvement in speed. We did, however, leave an option in settings for this stepping to be disabled completely, or for a custom value (e.g. every third layer) to be used.

File Name	Original (s)	Threaded (s)	Thread & Layer Optimised (s)
TOMO_IRNC001_DS.dcm	145.34	23.86	14.40
TOMO_IRChangAC001_DS.dcm	142.55	22.00	13.85

Table 2: Table of execution times for the uniformity calculations for provided DICOM images.

### 4.2 Code Structure Improvements

As visible in Appendix A, our code has three major Python files, `dicom_functions.py`, `uniformity_functions.py` and `gui.py`. Due to the complexity of the task and CustomTkinter naturally being quite verbose, this has resulted in some files of considerable length, for example `gui.py` at approximately 600 lines. While file length has no impact on the final result, this length combined with `gui.py` carrying out a large number of tasks made development difficult at times, as it was sometimes challenging to find (for example) the location of one function relative to the other. While we decided that overhauling the structure of the program late in development would not be a justified use of effort, if we had the opportunity to do this again we would pay more careful consideration to how we structured the project, for example a larger number of smaller, more modularised, files that carry out a much smaller number of tasks.

## 5 Acknowledgements

Section modified for anonymity in public copy

## 6 References

- [1] 3Blue1Brown. *But what is a convolution?* URL: <https://www.youtube.com/watch?v=KuXjwB4LzSA&t=1178s> (visited on 04/15/2024).
- [2] Robert Alley. *What Are Radioactive Tracers?* URL: <https://sciencing.com/radioactive-tracers-8330110.html> (visited on 04/15/2024).
- [3] National Institute of Biomedical Imaging and Engineering. *Nuclear Medicine*. URL: <https://www.nibib.nih.gov/science-education/science-topics/nuclear-medicine> (visited on 04/15/2024).
- [4] Wikipedia contributors. *Gamma camera – Wikipedia, The Free Encyclopedia*. URL: [https://en.wikipedia.org/w/index.php?title=Gamma\\_camera&oldid=1216441321](https://en.wikipedia.org/w/index.php?title=Gamma_camera&oldid=1216441321) (visited on 04/15/2024).
- [5] Wikipedia contributors. *Kernel (image processing) – Wikipedia, The Free Encyclopedia*. URL: [https://en.wikipedia.org/w/index.php?title=Kernel\\_\(image\\_processing\)&oldid=1180652863](https://en.wikipedia.org/w/index.php?title=Kernel_(image_processing)&oldid=1180652863) (visited on 04/15/2024).
- [6] DICOM. *6.2 Value Representation (VR)*. URL: [https://dicom.nema.org/medical/dicom/current/output/chtml/part05/sect\\_6.2.html#table\\_6.2-1](https://dicom.nema.org/medical/dicom/current/output/chtml/part05/sect_6.2.html#table_6.2-1) (visited on 04/15/2024).
- [7] DICOM. *8 Encoding of Pixel, Overlay and Waveform Data*. URL: [https://dicom.nema.org/medical/dicom/current/output/chtml/part05/chapter\\_8.html](https://dicom.nema.org/medical/dicom/current/output/chtml/part05/chapter_8.html) (visited on 04/17/2024).
- [8] DICOM. *PS3.6*. URL: <https://dicom.nema.org/medical/dicom/current/output/pdf/part06.pdf> (visited on 04/15/2024).
- [9] Sequoia Healthcare. *MiE Ecam Scintron Single Dual Head Gamma Camera System*. URL: <https://www.sqhpl.com/mie-ecam-scintron-single-dual-head-camma-camera.php> (visited on 04/15/2024).
- [10] LEAD Technologies Inc. *Overview: Basic DICOM File Structure*. URL: <https://www.leadtools.com/help/sdk/v21/dicom/api/overview-basic-dicom-file-structure.html> (visited on 04/15/2024).
- [11] LEAD Technologies Inc. *Overview: Data Element Structure*. URL: <https://www.leadtools.com/help/sdk/v21/dicom/api/overview-data-element-structure.html> (visited on 04/15/2024).
- [12] Michael K.O'Connor. *Quality Control of Scintillation Cameras (Planar and SPECT)*. URL: <https://www.aapm.org/meetings/99am/pdf/2741-51264.pdf> (visited on 04/19/2024).
- [13] Pypi. *customtkinter 5.2.2*. URL: <https://pypi.org/project/customtkinter/> (visited on 01/08/2024).
- [14] tbd. *Example of 2D Convolution*. URL: [http://www.songho.ca/dsp/convolution/convolution2d\\_example.html](http://www.songho.ca/dsp/convolution/convolution2d_example.html) (visited on 04/21/2024).
- [15] Cancer Research UK. *Radioactive scans for neuroendocrine tumours (NETs)*. URL: <https://www.cancerresearchuk.org/about-cancer/neuroendocrine-tumours-nets/getting-diagnosed/tests/radioactive-scans> (visited on 04/15/2024).



- 
- [16] Vaia. *Gamma Camera*. URL: <https://www.vaia.com/en-us/explanations/physics/medical-physics/gamma-camera/> (visited on 04/15/2024).

## A Code Appendix

<https://github.com/AlexToogood-Johnson/Gamma-Camera-Uniformity>

### A.1 dicom\_functions.py

```

1 # dicom_functions.py
2 """Functions relating to the manipulation of DICOM images"""
3
4 ##### IMPORTS #####
5
6 import numpy as np
7 import os
8 import json
9 from typing import NoReturn
10 from scipy.signal import convolve2d
11 import struct
12
13 ##### FUNCTIONS #####
14
15
16 def get_dicom_data(filepath: str) -> list:
17     """Reads the binary of a DICOM file"""
18
19     with open(filepath, "rb") as dicom_file:
20         data = dicom_file.read()
21     return data
22
23
24 def decode_value(vr, value):
25     """Decodes the value of a DICOM element based on its VR"""
26
27     try:
28         if vr in "AE AS CS DA DS DT LO LT PN SH ST TM UC UI UR UT UV":
29             return value.decode("ascii")
30         match vr:
31             case "AT": return value.hex()
32             case "FL": return struct.unpack(">f", value)[0]
33             case "FD": return struct.unpack(">d", value)[0]
34             case "IS": return int(value.decode("ascii"))
35             case "SL": return struct.unpack(">i", value)[0]
36             case "SS": return struct.unpack(">h", value)[0]
37             case "SV": return struct.unpack(">q", value)[0]
38             case "UL": return int.from_bytes(value, byteorder="little")
39             case "US": return int.from_bytes(value, byteorder="little")
40             case _: return value
41
42     except Exception:
43         return decode_value(vr, value[:-2])
44
45
46 def rearrange_tag(tag: str) -> str:
47     """When reading tags from binary, they are in the wrong order"""
48
49     tag = tag[2:4] + tag[0:2] + tag[6:] + tag[4:6]
50     return tag
51
52
53 def parse_binary(data) -> dict:

```

```

54     """Parses the raw binary from a DICOM file, and returns the metadata and raw
image data"""
55
56     parsed_data = {}
57     image_data = []
58     data = data[128:] # Remove preamble
59     data = data[4:] # Remove prefix
60
61     while len(data) > 0:
62         tag = rearrange_tag(data[:4].hex())
63         data = data[4:]
64         vr = data[:2].decode("ascii")
65         data = data[2:]
66         if vr in "AE AS AT CS DA DS DT FL FD IS LO LT PN SH SL ST SS TM UI UL US"
:
67             length = int.from_bytes(data[:2], byteorder="little")
68             data = data[2:]
69         else:
70             data = data[2:]
71             length = int.from_bytes(data[:4], byteorder="little")
72             data = data[4:]
73             value = decode_value(vr, data[:length])
74             data = data[length:]
75
76             if tag == "7fe00010":
77                 image_data = value
78                 data = []
79                 break
80             else:
81                 parsed_data[tag] = [vr, length, value]
82
83     elements = json.load(open("dicom_elements.json"))
84
85     for key, value in parsed_data.items():
86         if key in elements:
87             parsed_data[key].append(elements[key])
88         else:
89             parsed_data[key].append("UNKOWN")
90
91     return parsed_data, image_data
92
93
94 def decode_image_data(parsed_data: dict, image_data: list) -> np.ndarray:
95     """Converts the raw image data into a 3D Numpy array"""
96
97     bits_allocated = parsed_data["00280100"][2] // 8
98
99     image_data = [int.from_bytes(image_data[i:i + bits_allocated], byteorder="
little") for i in range(0, len(image_data), bits_allocated)]
100     rows = parsed_data["00280010"][2]
101     columns = parsed_data["00280011"][2]
102     planes = parsed_data["00280008"][2]
103     image_data = np.array(image_data).reshape((planes, rows, columns))
104
105     return image_data
106
107
108 def load_dicom_image(filepath: str) -> tuple:
109     """Returns the metadata and image data of a DICOM file"""
110
111     if not os.path.isfile(filepath):

```

```

112         raise FileNotFoundError(f"Error: {filepath} not found")
113     if not filepath.endswith(".dcm"):
114         raise ValueError(f"Error: {filepath} is not a DICOM file")
115
116     parsed_data, image_data = parse_binary(get_dicom_data(filepath))
117     image_data = decode_image_data(parsed_data, image_data)
118
119     return parsed_data, image_data
120
121
122 def read_config_file(key: str) -> list | NoReturn:
123     """Returns the contents of config.json"""
124
125     if not os.path.isfile("config.json"):
126         raise FileNotFoundError(f"Error: config.json not found, please run setup.
127 py")
128
129     with open("config.json", "r") as config_file:
130         data = json.load(config_file)
131     if key not in data:
132         raise KeyError(f"Error: key {key} not found in config.json")
133     return data[key]
134
135
136 def edit_config_file(key: str, value: list) -> NoReturn:
137     """Edits the contents of config.json"""
138
139     if not os.path.isfile("config.json"):
140         raise FileNotFoundError(f"Error: config.json not found, please run setup.
141 py")
142
143     with open("config.json", "r") as config_file:
144         data = json.load(config_file)
145     if key not in data:
146         raise KeyError(f"Error: key {key} not found in config.json")
147     data[key] = value
148     with open("config.json", "w") as config_file:
149         json.dump(data, config_file, indent=4)
150
151
152 def crop(array: np.ndarray, crop_size: int) -> np.ndarray | NoReturn:
153     """Returns the central n x n x n crop of the input array"""
154
155     if not len(array.shape) == 3:
156         raise ValueError(f"Expected 3D array. Got {len(array.shape)}D array
157 instead")
158
159     center = array.shape[0] // 2
160     start = center - crop_size // 2
161     end = start + crop_size
162     return array[start:end, start:end, start:end]
163
164
165 def apply_convolution(array: np.ndarray, convolution: list) -> np.ndarray |
166 NoReturn:
167     """Applies a smoothing convolution to a 2D Numpy array"""
168
169     if not len(array.shape) == 3:
170         raise ValueError(f"Expected 3D array. Got {len(array.shape)}D array
171 instead")
172     new_array = np.zeros_like(array)

```

```

168     for i in range(array.shape[0]):
169         new_array[i] = convolve2d(array[i], convolution, mode='same', boundary='
fill', fillvalue=0)
170     return new_array
171
172
173 def remove_image_edges(array: np.ndarray) -> np.ndarray | NoReturn:
174     """Removes the edges of a 3D Numpy array"""
175
176     if not len(array.shape) == 3:
177         raise ValueError(f"Expected 3D array. Got {len(array.shape)}D array
instead")
178
179     radius = array.shape[0] // 2
180     for layer in range(len(array)):
181         for row in range(len(array[layer])):
182             for column in range(len(array[layer][row])):
183                 if (radius - column) ** 2 + (radius - row) ** 2 > radius ** 2:
184                     array[layer][row][column] = 10000
185
186     return array
187
188
189 def max_pixel(array: list) -> int | NoReturn:
190     """Returns the maximum value in a list, ignoring None values"""
191
192     return max([value for value in array if value != 10000])
193
194
195 def min_pixel(array: list) -> int | NoReturn:
196     """Returns the minimum value in a list, ignoring None values"""
197
198     return min([value for value in array if value != 10000])

```

## A.2 uniformity\_functions.py

```

1 # uniformity_functions.py
2 """Functions relating to uniformity calculations"""
3
4 import numpy as np
5 from PIL import Image, ImageDraw
6 import concurrent.futures
7
8
9 class UniformityLayer:
10     def __init__(self, bin_data: list, circ_rad: int) -> None:
11         self.bin_data = np.array(bin_data, dtype=np.uint8)
12         self.circ_rad = circ_rad
13         self.cropped_data = None
14
15     def crop_to_circle(self) -> None:
16         image = Image.fromarray(self.bin_data)
17         if image.mode != 'RGBA':
18             image = image.convert('RGBA')
19         transparent_image = Image.new('RGBA', image.size)
20         transparent_image.paste(image, (0, 0))
21         mask = Image.new('L', image.size, 0)
22         draw = ImageDraw.Draw(mask)
23         center_x, center_y = image.width // 2, image.height // 2
24         draw.ellipse((center_x - self.circ_rad, center_y - self.circ_rad,

```

```

25         center_x + self.circ_rad, center_y + self.circ_rad), fill
=255)
26     transparent_image.putalpha(mask)
27     self.cropped_data = np.array(transparent_image)
28
29     def _calculate_uniformity(self, slices):
30         max_uniformity = 0
31         for s in slices:
32             if np.isnan(s).all():
33                 continue # Skip entirely transparent slices
34             valid_slice = s[~np.isnan(s)]
35             if valid_slice.size >= 2:
36                 max_brightness = np.nanmax(valid_slice)
37                 min_brightness = np.nanmin(valid_slice)
38                 if max_brightness + min_brightness != 0:
39                     uniformity = abs((max_brightness - min_brightness) / ((
max_brightness + min_brightness) / 2) * 100)
40                     max_uniformity = max(max_uniformity, uniformity)
41         return max_uniformity
42
43     def differential(self) -> float:
44         if self.cropped_data is None:
45             print("No data to process.")
46             return 0
47
48         alpha = self.cropped_data[:, :, 3] > 0
49         grayscale = 0.299 * self.cropped_data[:, :, 0] + \
50             0.587 * self.cropped_data[:, :, 1] + \
51             0.114 * self.cropped_data[:, :, 2]
52         grayscale[~alpha] = np.nan
53
54         horizontal_slices = [grayscale[y, x:x + 5] for y in range(grayscale.shape
[0]) for x in range(grayscale.shape[1] - 4)]
55         vertical_slices = [grayscale[y:y + 5, x] for x in range(grayscale.shape
[1]) for y in range(grayscale.shape[0] - 4)]
56
57         with concurrent.futures.ThreadPoolExecutor() as executor:
58             futures = [
59                 executor.submit(self._calculate_uniformity, horizontal_slices),
60                 executor.submit(self._calculate_uniformity, vertical_slices)
61             ]
62             results = [future.result() for future in concurrent.futures.
as_completed(futures)]
63
64         return max(results)
65
66     def integral(self) -> float:
67         if self.cropped_data is None:
68             print("No data to process.")
69             return 0
70
71         alpha = self.cropped_data[:, :, 3]
72         grayscale = 0.299 * self.cropped_data[:, :, 0] + \
73             0.587 * self.cropped_data[:, :, 1] + \
74             0.114 * self.cropped_data[:, :, 2]
75         grayscale[alpha == 0] = np.nan
76         valid_grayscale = grayscale[~np.isnan(grayscale)]
77
78         if valid_grayscale.size == 0:
79             return 0
80

```

```

81     max_brightness = np.nanmax(valid_grayscale)
82     min_brightness = np.nanmin(valid_grayscale)
83
84     integral_uniformity = abs((max_brightness - min_brightness) / ((
max_brightness + min_brightness) / 2) * 100)
85     return integral_uniformity

```

### A.3 gui.py

```

1  # gui.py
2  """The GUI for the gamma camera uniformity program
3  Copyright (c) 2024 Alex Toogood-Johnson, Lyall Stewart, Josh Scates, Leah
4  Wells, Zac Baker"""
5
6  ##### IMPORTS #####
7  import customtkinter as ctk
8  from tkinterdnd2 import TkinterDnD, DND_ALL
9  from PIL import Image, ImageTk, ImageDraw
10 import numpy as np
11 from dicom_functions import *
12 from uniformity_functions import UniformityLayer
13 from tkinter import filedialog
14 import os
15 import ast
16 from concurrent.futures import ProcessPoolExecutor
17 import time
18
19 ##### CONSTANTS #####
20
21 HEIGHT, WIDTH = 540, 960
22 FILEPATH = os.path.dirname(os.path.abspath(__file__))
23
24 ctk.set_default_color_theme(read_config_file("colour"))
25 ctk.set_appearance_mode(read_config_file("colour_theme"))
26
27 ##### CLASSES #####
28
29 class Tk(ctk.CTk, TkinterDnD.DnDWrapper):
30     """Sourced from: https://stackoverflow.com/questions/75526264/using-drag-and-drop-files-or-file-picker-with-customtkinter
31     Allows for drag and drop file opening within customtkinter"""
32
33     def __init__(self, *args, **kwargs) -> None:
34         super().__init__(*args, **kwargs)
35         self.TkdndVersion = TkinterDnD._require(self)
36
37
38 class Gui(Tk):
39     """The main GUI class for the Gamma Camera Uniformity program"""
40
41     def __init__(self, *args, **kwargs) -> None:
42         Tk.__init__(self, *args, **kwargs)
43
44         BUTTON_WIDTH = 150
45         TEXTBOX_WIDTH = 200
46         COMPONENT_HEIGHT = 30
47
48         self.bind('<Left>', lambda event: self.left_button_callback())
49

```

```

50     self.bind('<Right>', lambda event: self.right_button_callback())
51     self.bind('<Up>', lambda event: self.left_button_callback())
52     self.bind('<Down>', lambda event: self.right_button_callback())
53     self.bind('<o>', lambda event: self.open_button_callback())
54     self.bind('<s>', lambda event: self.apply_convolution())
55     self.bind('<BackSpace>', lambda event: self.revert_changes())
56
57     self.title("Gamma Camera Uniformity")
58     self.geometry(f"{WIDTH}x{HEIGHT}")
59     self.resizable(True, True)
60
61     self.photo_frame = ctk.CTkFrame(self, height=HEIGHT - 20, width=HEIGHT -
20)
62     self.photo_frame.place(x=10, y=10)
63     self.canvas = ctk.CTkCanvas(self.photo_frame, width=HEIGHT - 20, height=
HEIGHT - 20)
64     self.canvas.place(x=0, y=0)
65     self.canvas_image = None
66     self.tk_image = None
67
68     self.dnd_box = ctk.CTkEntry(self.photo_frame, width=HEIGHT - 20, height=
HEIGHT - 20, state="readonly")
69     self.dnd_box.place(x=0, y=0)
70     self.dnd_box.drop_target_register(DND_ALL)
71     self.dnd_box.dnd_bind("<Drop>", self.get_path)
72
73     self.tab_view = ctk.CTkTabview(self, width=WIDTH - HEIGHT - 10, height=
HEIGHT - 10)
74     self.tab_view.place(x=HEIGHT, y=0)
75     self.tab_1 = self.tab_view.add("    DICOM Operations    ")
76     self.tab_2 = self.tab_view.add("    Uniformity Calculations    ")
77     self.tab_3 = self.tab_view.add("    Settings    ")
78
79     self.frame_tab_1 = ctk.CTkFrame(self.tab_1, width=WIDTH - HEIGHT - 20)
80     self.frame_tab_1.pack(side='top', pady=40)
81     self.combobox_values = [str(i) for i in range(1, 40)]
82     self.left_button = ctk.CTkButton(self.frame_tab_1, text="Left", width=70,
height=10, command=self.left_button_callback, font=("", 20), state="disabled"
)
83     self.select_layer = ctk.CTkComboBox(self.frame_tab_1, width=190, height
=10, state="disabled", font=("", 20), values=self.combobox_values, command=
self.combobox_callback)
84     self.right_button = ctk.CTkButton(self.frame_tab_1, text="Right", width
=70, height=10, command=self.right_button_callback, font=("", 20), state="
disabled")
85     self.left_button.pack(side='left', padx=10)
86     self.select_layer.pack(side='left', padx=10)
87     self.right_button.pack(side='left', padx=10)
88
89     self.frame_tab_2 = ctk.CTkFrame(self.tab_1, width=WIDTH - HEIGHT - 20)
90     self.frame_tab_2.pack()
91     self.open_button = ctk.CTkButton(self.frame_tab_2, text="    Open    ",
width=30, height=20, font=("", 20), command=self.open_button_callback)
92     self.close_button = ctk.CTkButton(self.frame_tab_2, text="    Close    ",
width=30, height=20, font=("", 20), state="disabled", command=self.
close_button_callback)
93     self.save_button = ctk.CTkButton(self.frame_tab_2, text="    Save Layer    ",
width=30, height=20, font=("", 20), state="disabled", command=self.
save_button_callback)
94     self.open_button.pack(side='left', padx=10)
95     self.save_button.pack(side='left', padx=10)

```



```

96         self.close_button.pack(side='left', padx=10)
97         self.open_button.configure(state='normal')
98
99         self.frame_tab_4 = ctk.CTkFrame(self.tab_1, width=WIDTH - HEIGHT - 20)
100         self.frame_tab_4.pack(pady=40)
101         self.apply_convolution_button = ctk.CTkButton(self.frame_tab_4, text="
Smooth ", width=30, height=20, font=("", 20), state="disabled", command=self.
apply_convolution)
102         self.reduce_image_button = ctk.CTkButton(self.frame_tab_4, text="
Reduce ", width=30, height=20, font=("", 20), state="disabled", command=self.
.reduce_image)
103         self.revert_changes_button = ctk.CTkButton(self.frame_tab_4, text="
Revert ", width=30, height=20, font=("", 20), state="disabled", command=self.
.revert_changes)
104         self.apply_convolution_button.pack(side='left', padx=10)
105         self.reduce_image_button.pack(side='left', padx=10)
106         self.revert_changes_button.pack(side='left', padx=10)
107
108         self.frame_tab_5 = ctk.CTkFrame(self.tab_1, width=WIDTH - HEIGHT - 20,
height=190)
109         self.frame_tab_5.pack(side='bottom', pady=20)
110         self.logo = Image.open(os.path.join(FILEPATH, "logo.png"))
111         self.logo = self.logo.resize((WIDTH - HEIGHT - 20, 190), Image.LANCZOS)
112         self.tk_logo = ImageTk.PhotoImage(self.logo)
113         self.logo_label = ctk.CTkLabel(self.frame_tab_5, image=self.tk_logo, text
="")
114         self.logo_label.pack()
115
116         self.frame_tab_6 = ctk.CTkFrame(self.tab_3, width=WIDTH - HEIGHT - 20)
117         self.frame_tab_6.pack(pady=10)
118         self.about_button = ctk.CTkButton(self.frame_tab_6, text=" About ",
width=30, height=20, font=("", 20), command=lambda: AboutGUI().mainloop())
119         self.help_button = ctk.CTkButton(self.frame_tab_6, text=" Help ",
width=30, height=20, font=("", 20), command=lambda: HelpGUI().mainloop())
120         self.licence_button = ctk.CTkButton(self.frame_tab_6, text=" Licence
", width=30, height=20, font=("", 20), command=lambda: LicenceGUI().mainloop()
)
121         self.about_button.pack(side='left', padx=10)
122         self.help_button.pack(side='left', padx=10)
123         self.licence_button.pack(side='left', padx=10)
124
125         self.frame_tab_7 = ctk.CTkFrame(self.tab_3, width=WIDTH - HEIGHT - 20)
126         self.frame_tab_7.pack(pady=(10, 2), anchor='w')
127         self.change_saving_directory_button = ctk.CTkButton(self.frame_tab_7,
text="Save To", width=30, height=30, font=("", 15), command=self.
save_directory_button_callback)
128         dirname = read_config_file("default_file_saving_directory")
129         self.saving_directory_name = ctk.CTkLabel(self.frame_tab_7, text=self.
normalize_directory(dirname), font=("courier", 15))
130         self.change_saving_directory_button.pack(side='left', padx=10)
131         self.saving_directory_name.pack(side='left', padx=10, pady=5)
132
133         self.frame_tab_8 = ctk.CTkFrame(self.tab_3, width=WIDTH - HEIGHT - 20)
134         self.frame_tab_8.pack(pady=(2, 10), anchor='w')
135         self.change_opening_directory_button = ctk.CTkButton(self.frame_tab_8,
text="Open At", width=30, height=30, font=("", 15), command=self.
open_directory_button_callback)
136         dirname = read_config_file("default_file_opening_directory")
137         self.opening_directory_name = ctk.CTkLabel(self.frame_tab_8, text=self.
normalize_directory(dirname), font=("courier", 15))
138         self.change_opening_directory_button.pack(side='left', padx=10)

```

```

139         self.opening_directory_name.pack(side='left', padx=10, pady=10)
140
141         self.frame_tab_9 = ctk.CTkFrame(self.tab_3, width=WIDTH - HEIGHT - 20)
142         self.frame_tab_9.pack(pady=10, anchor='w')
143         self.theme_menu = ctk.CTkOptionMenu(self.frame_tab_9, font=("", 15),
width=BUTTON_WIDTH, command=lambda event: self.theme_changed(event), values=["
Light", "Dark", "System"])
144         self.theme_menu.pack(side='left', padx=10)
145         self.theme_menu.set(read_config_file("colour_theme"))
146         self.colour_menu = ctk.CTkOptionMenu(self.frame_tab_9, font=("", 15),
width=BUTTON_WIDTH, command=lambda event: self.colour_changed(event), values=["
blue", "green", "dark-blue"])
147         self.colour_menu.pack(side='left', padx=10)
148         self.colour_menu.set(read_config_file("colour"))
149
150         self.frame_tab_10 = ctk.CTkFrame(self.tab_3, width=WIDTH-HEIGHT-20)
151         self.frame_tab_10.pack(pady=5, anchor='w')
152
153         self.change_default_crop_size_button = ctk.CTkButton(self.frame_tab_10,
text="Update Crop Size", width=BUTTON_WIDTH, height=COMPONENT_HEIGHT, font=("",
, 15), command=self.change_crop_size_callback)
154         self.change_default_crop_size_button.pack(side='left', padx=10)
155         self.crop_size_textbox = ctk.CTkTextbox(self.frame_tab_10, width=
TEXTBOX_WIDTH, height=1, font=("", 15))
156         self.crop_size_textbox.insert("1.0", str(read_config_file("crop_amount"))
)
157         self.crop_size_textbox.pack(side='left', padx=10)
158
159         self.frame_tab_11 = ctk.CTkFrame(self.tab_3, width=WIDTH - HEIGHT - 20)
160         self.frame_tab_11.pack(pady=(10, 0), anchor='w')
161         self.change_fov_radius_button = ctk.CTkButton(self.frame_tab_11, text="
Update FoV Radius", width=BUTTON_WIDTH, height=COMPONENT_HEIGHT, font=("", 15)
, command=self.change_fov_radius_callback)
162         self.change_fov_radius_button.pack(side='left', padx=10)
163         self.fov_radius_textbox = ctk.CTkTextbox(self.frame_tab_11, width=
TEXTBOX_WIDTH, height=1, font=("", 15))
164         self.fov_radius_textbox.insert("1.0", str(read_config_file("fov_radius"))
)
165
166         self.fov_radius_textbox.pack(side='left', padx=10, pady=10)
167
168         self.frame_tab_14 = ctk.CTkFrame(self.tab_3, width=WIDTH-HEIGHT-20)
169         self.frame_tab_14.pack(pady=(0, 10), anchor='w')
170         self.draw_crop_size_button = ctk.CTkButton(self.frame_tab_14, text="Draw
FoV Radius", width=BUTTON_WIDTH + TEXTBOX_WIDTH + 20, height=COMPONENT_HEIGHT,
font=("", 15), command=self.draw_fov_radius_callback)
171         self.draw_crop_size_button.pack(side='left', padx=10, pady=10)
172         self.draw_crop_size_button.configure(state="disabled")
173
174         self.frame_tab_13 = ctk.CTkFrame(self.tab_3, width=WIDTH-HEIGHT-20)
175         self.frame_tab_13.pack(anchor='w')
176
177         self.change_step_button = ctk.CTkButton(self.frame_tab_13, text="Update
Step", width=BUTTON_WIDTH, height=COMPONENT_HEIGHT, font=("", 15), command=
self.change_step)
178         self.change_step_button.pack(side='left', padx=10)
179         self.step_textbox = ctk.CTkTextbox(self.frame_tab_13, width=TEXTBOX_WIDTH
, height=1, font=("", 15))
180         self.step_textbox.insert("1.0", str(read_config_file("step")))
181         self.step_textbox.pack(side='left', padx=10)
182

```

```

183     self.frame_tab_15 = ctk.CTkFrame(self.tab_3, width=WIDTH - HEIGHT - 20)
184     self.frame_tab_15.pack(pady=10, anchor='w')
185     self.reset_button = ctk.CTkButton(self.frame_tab_15, text="Reset to
recommended settings", width=BUTTON_WIDTH + TEXTBOX_WIDTH + 20, height=30,
font=("", 15), command=self.reset_settings)
186     self.reset_button.pack(side='left', padx=10, pady=10)
187
188     self.differential_uniformity_frame = ctk.CTkFrame(self.tab_2, width=WIDTH
- HEIGHT - 20)
189     self.differential_uniformity_frame.pack(side='top', pady=40, anchor='w')
190
191     self.get_differential_uniformity_button = ctk.CTkButton(self.
differential_uniformity_frame, text="Get Uniformity", width=30, height=20,
font=("", 15), command=self.uniformity_callback)
192     self.get_differential_uniformity_button.pack(side='left', padx=10, anchor
='w')
193     self.differential_uniformity_label = ctk.CTkLabel(self.
differential_uniformity_frame, text="Differential Uniformity: ", font=("", 15)
)
194     self.differential_uniformity_label.pack(side='left', padx=10)
195     self.differential_text_output = ctk.CTkTextbox(self.tab_2, width=30,
height=10, font=("", 15))
196     self.differential_text_output.pack(padx=10, pady=10, expand=True, fill='
both')
197
198     self.integral_uniformity_frame = ctk.CTkFrame(self.tab_2, width=WIDTH -
HEIGHT - 20)
199     self.integral_uniformity_frame.pack(pady=40, anchor='w')
200     self.get_integral_uniformity_button = ctk.CTkButton(self.
integral_uniformity_frame, text="Get Uniformity", width=30, height=20, font=("",
15), command=self.uniformity_callback)
201     self.get_integral_uniformity_button.pack(side='left', padx=10, anchor='w'
)
202     self.integral_uniformity_label = ctk.CTkLabel(self.
integral_uniformity_frame, text="Integral Uniformity: ", font=("", 15))
203     self.integral_uniformity_label.pack(side='left', padx=10)
204     self.integral_text_output = ctk.CTkTextbox(self.tab_2, width=30, height
=10, font=("", 15))
205     self.integral_text_output.pack(padx=10, pady=10, expand=True, fill='both'
)
206
207     self.get_differential_uniformity_button.configure(state="disabled")
208     self.get_integral_uniformity_button.configure(state="disabled")
209
210     def change_step(self) -> None:
211         step_size = self.step_textbox.get("1.0", "end-1c")
212         if step_size.isdigit() and int(step_size) > 0 and int(step_size) < 40:
213             edit_config_file("crop", int(step_size))
214         else:
215             self.step_size.delete("1.0", "end")
216             self.step_size.insert("1.0", str(read_config_file("step")))
217
218     def reset_settings(self) -> None:
219         self.crop_size_textbox.delete("1.0", "end")
220         self.crop_size_textbox.insert("1.0", "40")
221         self.fov_radius_textbox.delete("1.0", "end")
222         self.fov_radius_textbox.insert("1.0", "80")
223         self.step_textbox.delete("1.0", "end")
224         self.step_textbox.insert("1.0", "2")
225         self.change_crop_size_callback()
226

```

```

227     def uniformity_callback(self) -> None:
228         self.cropped_dicom_image = apply_convolution(self.current_dicom_image,
229 read_config_file("convolution"))
230         self.cropped_dicom_image = crop(self.current_dicom_image,
231 read_config_file("crop_amount"))
232
233         self.differential_text_output.insert("1.0", "Loading...")
234         self.integral_text_output.insert("1.0", "Loading...")
235
236         differential_vals = []
237         integral_vals = []
238
239         with ProcessPoolExecutor() as executor:
240             future_to_type = {}
241             for layer_index in range(0, read_config_file("crop_amount") - 1, 2):
242                 dicom_slice = self.current_dicom_image[layer_index]
243                 normalized_slice = ((dicom_slice - np.min(dicom_slice)) / (np.max
244 (dicom_slice) - np.min(dicom_slice)) * 255).astype(np.uint8)
245
246                 image = Image.fromarray(normalized_slice)
247                 image = image.convert("RGB")
248                 image = image.resize((HEIGHT - 20, HEIGHT - 20), Image.LANCZOS)
249
250                 layer = UniformityLayer(np.array(image), read_config_file("
251 fov_radius"))
252                 layer.crop_to_circle()
253
254                 future_to_type[executor.submit(layer.differential)] = '
255 differential'
256                 future_to_type[executor.submit(layer.integral)] = 'integral'
257
258                 start_time = time.time()
259                 for future, type_label in future_to_type.items():
260                     result = future.result()
261                     if type_label == 'differential':
262                         differential_vals.append(result)
263                     else:
264                         integral_vals.append(result)
265                 end_time = time.time()
266
267                 if differential_vals:
268                     self.display_uniformity(integral_vals, differential_vals)
269                 else:
270                     print("No differential values to display.")
271
272                 print(f"Total Execution Time: {end_time - start_time} seconds")
273
274     def display_uniformity(self, integral: int, differential: int) -> None:
275         self.differential_uniformity_results = differential
276         self.differential_uniformity_label.configure(text=f" Differential
277 Uniformity: {round(max(self.differential_uniformity_results), 2)}%")
278         fov = read_config_file("fov_radius")
279         convolution = read_config_file("convolution")
280         textbox_text = f"Differential Uniformity Results:\nFoV Radius {fov} px \
281 nConvolution: {convolution}\n\n"
282         for i, result in enumerate(self.differential_uniformity_results):
283             textbox_text += f"Layer {(i+1) * read_config_file('step')}: {result}\
284 n"
285
286         self.differential_text_output.delete("1.0", "end")
287         self.differential_text_output.insert("1.0", textbox_text)
288
289

```

```

280     self.integral_uniformity_results = integral
281     self.integral_uniformity_label.configure(text=f"Integral Uniformity: {
round(max(self.integral_uniformity_results), 2)}%")
282     fov = read_config_file("fov_radius")
283     convolution = read_config_file("convolution")
284     textbox_text = f"Integral Uniformity Results:\nFoV Radius {fov} px \
nConvolution: {convolution}\n\n"
285     for i, result in enumerate(self.integral_uniformity_results):
286         textbox_text += f"Layer {(i+1) * read_config_file('step')}: {result}\
n"
287     self.integral_text_output.delete("1.0", "end")
288     self.integral_text_output.insert("1.0", textbox_text)
289
290     def change_crop_size_callback(self) -> None:
291         """Changes the crop size stored in the config file"""
292
293         crop_size = self.crop_size_textbox.get("1.0", "end-1c")
294         if crop_size.isdigit() and int(crop_size) > 0 and int(crop_size) < 100:
295             edit_config_file("crop_amount", int(crop_size))
296         else:
297             self.crop_size_textbox.delete("1.0", "end")
298             self.crop_size_textbox.insert("1.0", str(read_config_file("
crop_amount"))))
299
300     def change_fov_radius_callback(self) -> None:
301         """Changes the field of view radius stored in the config file"""
302         fov_rad = self.fov_radius_textbox.get("1.0", "end-1c")
303         if fov_rad.isdigit() and int(fov_rad) > 0 and int(fov_rad) < 100:
304             edit_config_file("fov_radius", int(fov_rad))
305         else:
306             self.fov_radius_textbox.delete("1.0", "end")
307             self.fov_radius_textbox.insert("1.0", str(read_config_file("
fov_radius"))))
308
309     def draw_fov_radius_callback(self) -> None:
310         """Draws the field of view radius on the image"""
311         self.display_image(int(self.select_layer.get()), int(self.
fov_radius_textbox.get("1.0", "end")))
312
313     def colour_changed(self, event) -> None:
314         """Changes the colour of the GUI"""
315         colour = self.colour_menu.get()
316         ctk.set_default_color_theme(colour)
317         edit_config_file("colour", colour)
318
319     def theme_changed(self, event) -> None:
320         """Changes the colour theme of the GUI"""
321         theme = self.theme_menu.get()
322         ctk.set_appearance_mode(theme)
323         edit_config_file("colour_theme", theme)
324
325     def save_directory_button_callback(self) -> None:
326         """Changes the directory stored in the config file for saving files"""
327         directory = filedialog.askdirectory()
328         if directory:
329             if os.path.isdir(directory):
330                 self.saving_directory_name.configure(text=self.
normalize_directory(directory))
331                 edit_config_file("default_file_saving_directory", directory)
332
333     def normalize_directory(self, directory: str) -> str:

```

```

334     """Makes all directory names 30 character long for display purposes"""
335     if len(directory) >= 25:
336         directory = "... " + directory[-25:] + " "
337     else:
338         directory = " " + directory + " " + " " * (28 - len(directory))
339     return directory
340
341     def open_directory_button_callback(self) -> None:
342         """Changes the directory stored in the config file for opening files"""
343         directory = filedialog.askdirectory()
344         if directory:
345             if os.path.isdir(directory):
346                 self.opening_directory_name.configure(text=self.
normalize_directory(directory))
347                 edit_config_file("default_file_opening_directory", directory)
348
349     def reduce_image(self) -> None:
350         """Reduces the image by the amount specified in the config file, and
displays in the GUI"""
351         if self.reduce_image_button.cget("state") != 'normal': return
352         self.current_dicom_image = crop(self.current_dicom_image,
read_config_file("crop_amount"))
353         self.reduce_image_button.configure(state="disabled")
354         self.apply_convolution_button.configure(state="disabled")
355         self.revert_changes_button.configure(state="normal")
356         self.select_layer.set("1")
357         self.left_button.configure(state="disabled")
358         self.combobox_values = [str(i) for i in range(1, read_config_file("
crop_amount") - 1)]
359         self.display_image(int(self.select_layer.get()))
360
361     def apply_convolution(self) -> None:
362         """Applies a convolution to the image, and displays in the GUI"""
363         if self.apply_convolution_button.cget("state") != 'normal': return
364         self.current_dicom_image = apply_convolution(self.current_dicom_image,
read_config_file("convolution"))
365         self.display_image(int(self.select_layer.get()))
366         self.apply_convolution_button.configure(state="disabled")
367         self.reduce_image_button.configure(state="normal")
368         self.revert_changes_button.configure(state="normal")
369
370     def revert_changes(self) -> None:
371         """Reverts the visual changes made to the image - convolution and/or
cropping"""
372         if self.revert_changes_button.cget("state") != 'normal': return
373         self.current_dicom_image = self.original_dicom_image
374         self.display_image(int(self.select_layer.get()))
375         self.reduce_image_button.configure(state="disabled")
376         self.apply_convolution_button.configure(state="normal")
377         self.revert_changes_button.configure(state="disabled")
378
379     def combobox_callback(self, choice: str) -> None:
380         """Callback for selection of layer number from combobox"""
381         if int(choice) == 1:
382             self.left_button.configure(state="disabled")
383             self.right_button.configure(state="normal")
384         elif choice == self.combobox_values[-1]:
385             self.left_button.configure(state="normal")
386             self.right_button.configure(state="disabled")
387         else:
388             self.left_button.configure(state="normal")

```



```

389         self.right_button.configure(state="normal")
390         self.display_image(int(choice))
391
392     def left_button_callback(self) -> None:
393         """Callback for left button press on layer selection"""
394         if self.left_button.cget("state") != 'normal': return
395         self.select_layer.set(int(self.select_layer.get()) - 1)
396         if int(self.select_layer.get()) == 1:
397             self.left_button.configure(state="disabled")
398         if self.select_layer.get() != self.combobox_values[-1]:
399             self.right_button.configure(state="normal")
400         self.display_image(int(self.select_layer.get()))
401
402     def right_button_callback(self):
403         """Callback for the right button press on layer selection"""
404         if self.right_button.cget("state") != 'normal': return
405         self.select_layer.set(int(self.select_layer.get()) + 1)
406         if int(self.select_layer.get()) > 1:
407             self.left_button.configure(state="normal")
408         if self.select_layer.get() == self.combobox_values[-1]:
409             self.right_button.configure(state="disabled")
410         self.display_image(int(self.select_layer.get()))
411
412     def open_button_callback(self) -> None:
413         """Callback for the open button press, opens a file dialog to select a
414         DICOM file to open"""
415         if self.open_button.cget("state") != 'normal': return
416         default = read_config_file("default_file_opening_directory")
417         if os.path.isdir(default):
418             file_path = filedialog.askopenfilename(filetypes=[("DICOM Files", "*.
419             dcm")], initialdir=default)
420         else:
421             file_path = filedialog.askopenfilename(filetypes=[("DICOM Files", "*.
422             dcm")], initialdir=os.getcwd())
423         if file_path:
424             self.title(f"Gamma Camera Uniformity - {file_path.split('/')[-1]}")
425             self.fit_dicom_image(file_path)
426
427     def close_button_callback(self) -> None:
428         """Callback for the close button press, resets the GUI to its initial
429         state"""
430         self.current_dicom_image = None
431         self.original_dicom_image = None
432         self.combobox_values = []
433         self.select_layer.set("")
434         self.left_button.configure(state="disabled")
435         self.right_button.configure(state="disabled")
436         self.select_layer.configure(state="disabled")
437         self.open_button.configure(state="normal")
438         self.save_button.configure(state="disabled")
439         self.close_button.configure(state="disabled")
440         self.reduce_image_button.configure(state="disabled")
441         self.apply_convolution_button.configure(state="disabled")
442         self.revert_changes_button.configure(state="disabled")
443         self.get_integral_uniformity_button.configure(state="disabled")
444         self.get_differential_uniformity_button.configure(state="disabled")
445         self.differential_text_output.delete("1.0", "end")
446         self.integral_text_output.delete("1.0", "end")
447         self.differential_uniformity_label.configure(text="Uniformity: ")
448         self.integral_uniformity_label.configure(text="Uniformity: ")
449         self.title("Gamma Camera Uniformity")

```

```

446         self.draw_crop_size_button.configure(state="disabled")
447         self.canvas.delete(self.canvas_image)
448         self.dnd_box.lift()
449
450     def save_button_callback(self) -> None:
451         """Callback for the save button press, saves the currently viewed layer
452         as a PNG file"""
453         dicom_slice = self.current_dicom_image[int(self.select_layer.get()) - 1]
454         normalized_slice = ((dicom_slice - np.min(dicom_slice)) / (np.max(
455             dicom_slice) - np.min(dicom_slice)) * 255).astype(np.uint8)
456         default = read_config_file("default_file_saving_directory")
457         if os.path.isdir(default):
458             file_path = filedialog.asksaveasfilename(defaultextension=".png",
459                 filetypes=[("PNG Files", "*.png")], initialdir=default)
460         else:
461             file_path = filedialog.asksaveasfilename(defaultextension=".png",
462                 filetypes=[("PNG Files", "*.png")], initialdir=os.getcwd())
463         if file_path:
464             Image.fromarray(normalized_slice).save(file_path)
465
466     def get_path(self, event) -> None:
467         image_path = event.data.replace("{", "").replace("}", "")
468         if image_path.endswith(".dcm"):
469             self.fit_dicom_image(image_path)
470
471     def fit_dicom_image(self, image_path) -> None:
472         self.current_dicom_image = load_dicom_image(image_path)[1]
473         if not len(self.current_dicom_image.shape) == 3:
474             self.current_dicom_image = None
475             return # Rejects 2D images
476         self.original_dicom_image = self.current_dicom_image
477         self.combobox_values = [str(i) for i in range(1, self.current_dicom_image
478             .shape[0] + 1)]
479         self.right_button.configure(state="normal")
480         self.select_layer.configure(state="normal")
481         self.select_layer.set("1")
482         self.open_button.configure(state="disabled")
483         self.save_button.configure(state="normal")
484         self.close_button.configure(state="normal")
485         self.apply_convolution_button.configure(state="normal")
486         self.get_integral_uniformity_button.configure(state="normal")
487         self.get_differential_uniformity_button.configure(state="normal")
488         self.draw_crop_size_button.configure(state="normal")
489         self.display_image(1)
490
491     def display_image(self, layer: int, fov_rad: int = 0) -> None:
492         dicom_slice = self.current_dicom_image[layer - 1]
493         normalized_slice = ((dicom_slice - np.min(dicom_slice)) / (np.max(
494             dicom_slice) - np.min(dicom_slice)) * 255).astype(np.uint8)
495
496         image = Image.fromarray(normalized_slice)
497         image = image.convert("RGB")
498         image = image.resize((HEIGHT - 20, HEIGHT - 20), Image.LANCZOS)
499         if fov_rad > 0: # Providing a FoV of 0 indicates that the user does not
500             want a FoV drawn
501             draw = ImageDraw.Draw(image)
502             center_x, center_y = (HEIGHT - 20) // 2, (HEIGHT - 20) // 2
503             top_left = (center_x - fov_rad, center_y - fov_rad)
504             bottom_right = (center_x + fov_rad, center_y + fov_rad)
505             draw.ellipse([top_left, bottom_right], outline="red", width=2)

```



```

500         self.tk_image = ImageTk.PhotoImage(image)
501         if self.canvas_image:
502             self.canvas.delete(self.canvas_image)
503         self.canvas_image = self.canvas.create_image(-1, -1, anchor="nw", image=
self.tk_image, tags="image")
504         self.canvas.lift(self.canvas_image)
505         self.dnd_box.lower()
506
507
508 class LicenceGUI(Tk):
509     def __init__(self, *args, **kwargs) -> None:
510         Tk.__init__(self, *args, **kwargs)
511
512         self.title("Licence")
513         self.geometry("500x300")
514         self.resizable(False, False)
515
516         self.licence_text = ctk.CTkTextbox(self, width=480, height=280)
517         self.licence_text.place(x=10, y=10)
518         self.licence_text.insert("1.0", open(os.path.join(FILEPATH, "LICENCE")).
read())
519         self.licence_text.configure(state="disabled")
520
521
522 class AboutGUI(Tk):
523     def __init__(self, *args, **kwargs) -> None:
524         Tk.__init__(self, *args, **kwargs)
525         self.title("About")
526         self.geometry("500x300")
527         self.resizable(False, False)
528         text = """
529 Gamma Camera Uniformity Program
530 Developed at Exeter Maths School as part of the Exeter Maths Certificate
531 in 2024.
532 The purpose of this software is to perform calculations to determine the
533 uniformity values of DICOM SPECT images taken from a gamma camera.
534
535 For more information, and to view the source code, visit the GitHub page:
536 https://github.com/ExeMS/NHS-EMC2024
537
538 This software is distributed under the MIT License.
539 Copyright (c) 2024 Alex Toogood-Johnson, Lyall Stewart, Josh Scates, Leah Wells,
Zac Baker.
540
541         """
542         self.about_text = ctk.CTkTextbox(self, width=480, height=280)
543         self.about_text.place(x=10, y=10)
544         self.about_text.insert("1.0", text)
545         self.about_text.configure(state="disabled")
546
547
548 class HelpGUI(Tk):
549     def __init__(self, *args, **kwargs) -> None:
550         Tk.__init__(self, *args, **kwargs)
551         self.title("Help")
552         self.geometry("500x300")
553         self.resizable(False, False)
554         text = """
555 This section provides some information on how to use this program, although for
further details please contact the authors. To start, open a DICOM file by
selecting it from the file directory, which you can access by typing 'o' on

```

your keyboard or by clicking the 'open' button. You can also directly open a file by holding it over the left half of the program, which functions as a drag and drop box as long as there isn't currently a file open.

The chosen file must be a .dcm file, and must be a SPECT image in 3 dimensions in order to be opened by the program. Once an image has been opened, you can close the image by clicking the 'close' button. In order to view each layer of the DICOM image, you can either click the left, right, up or down buttons on your keyboard, select the left or right buttons on the gui, or manually select a layer from the combobox. Please note that the first and last 20-40 layers in a DICOM uniformity image may be blank or incomplete.

The 'reduce' button, which you can either click, or press 'r' on your keyboard, fixes this by removing the first n and last n layers of the image. The default and recommended value is 44, although this can be changed in the settings tab.

The smooth button, which you can either click or press 's' on your keyboard, applies a 9 point weighted convolution in order to smooth each layer of the DICOM image. The convolution used is  $\begin{bmatrix} 1 & 2 & 1 \\ 2 & 4 & 2 \\ 1 & 2 & 1 \end{bmatrix}$  as specified in the NEMA handbook. We strongly recommend not changing this, however if necessary than this can be manually changed inside 'config.json'.

The revert button, which you can either click, or press backspace on your keyboard, visually reverts the smoothing and / or cropping changes to the DICOM image. The save layer button saves the currently selected layer of the DICOM image as a PNG file. The directory at which this opens, as well as the directory which is opened when you want to open a DICOM file, can be specified in the settings tab.

In the middle tab of the program, there are options to calculate integral and differential uniformity, although please note that pressing one of these buttons also triggers the other. Due to a high number of calculations in order to find the uniformities, this may take up to a minute, during which time the program may not respond. In the labels you will be able to see the overall differential or integral uniformity, whereas in the larger text box there will be details about each layer in the DICOM image. In the third tab, which contains the settings for the program, there are buttons linking to the about, help and licence sections. Add a bit more stuff here.

```
self.about_text = ctk.CTkTextbox(self, width=480, height=280)
self.about_text.place(x=10, y=10)
self.about_text.insert("1.0", text)
self.about_text.configure(state="disabled")
```

```
##### MAIN #####
```

```
if __name__ == "__main__":
    app = Gui()
    app.mainloop()
```

## A.4 config.json

```
{
  "default_file_opening_directory": "C:/Users/Alex Toogood-Johnson/Documents",
  "default_file_saving_directory": "C:/Users/Alex Toogood-Johnson/Downloads",
  "colour_theme": "System",
  "colour": "blue",
  "convolution": [
```

```
7      [
8          1,
9          2,
10         1
11     ],
12     [
13         2,
14         4,
15         2
16     ],
17     [
18         1,
19         2,
20         1
21     ]
22 ],
23 "crop_amount": 44,
24 "fov_radius": 80,
25 "step": 2
26 }
```

### A.5 dicom\_elements.json

A small sample of this file has been included to demonstrate its format, however due to its length of approximately 5000 lines, it has been heavily truncated.

```
1 {
2     "00020000": "File Meta Information Group Length",
3     "00020001": "File Meta Information Version",
4     "00020002": "Media Storage SOP Class UID",
5     "00020003": "Media Storage SOP Instance UID",
6     "00020010": "Transfer Syntax UID",
7     "00020012": "Implementation Class UID",
8     "00020013": "Implementation Version Name",
9     "00020016": "Source Application Entity Title",
10    ...
11 }
```